

Product datasheet for **AM08207FC-N**

Mouse IgG1 Isotype Control

Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody is specific to T-2 Mycotoxin.
Buffer:	1.0 mL of PBS containing 0.09% Sodium Azide as preservative.
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)