

## **Product datasheet for AM08207PU-N**

## **Mouse IgG1 Isotype Control**

## **Product data:**

**Product Type:** Isotype Controls

Applications: E, FC
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Specificity:** This antibody is specific to T-2 Mycotoxin.

**Buffer:** Borate buffered saline, pH 8.2.

No preservatives or amine-containing buffer salts added.

Concentration: 1 mg/ml

Conjugation: Unconjugated

**Storage:** Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing. Shelf life: one year from despatch.



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com