

Product datasheet for **AM08208AC5-N**

Mouse IgG2a Isotype Control

Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Buffer:	1.0 mL of PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.
Conjugation:	APC-Cy5.5
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)