

Product datasheet for **AM08208BT-N**

Mouse IgG2a Isotype Control

Product data:

Product Type: Isotype Controls

Applications: E, FC

Host: Mouse

Isotype: IgG2a

Clonality: Monoclonal

Buffer: PBS containing 0.09% Sodium Azide as preservative.

Concentration: 0.5 mg/ml

Conjugation: Biotin

Storage: Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)