

Product datasheet for **AM08208PU-N**

Mouse IgG2a Isotype Control

Product data:

Product Type: Isotype Controls

Applications: E, FC

Host: Mouse

Isotype: IgG2a

Clonality: Monoclonal

Buffer: 100 mM Borate buffered saline, pH 8.2.
No preservatives or amine-containing buffer salts added.

Concentration: 1 mg/ml

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.
Avoid repeated freezing and thawing.
Shelf life: one year from despatch.



[View online »](#)