

Product datasheet for **AM08209LE-N**

Mouse IgG2b Isotype Control

Product data:

Product Type:	Isotype Controls
Applications:	E, FC, ICC/IF
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Specificity:	This antibody is specific for Chicken IgA.
Buffer:	PBS without preservatives.
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)