

## Product datasheet for **AM08210FC-N**

### Mouse IgG3 Isotype Control

#### Product data:

**Product Type:** Isotype Controls

**Applications:** FC

**Host:** Mouse

**Isotype:** IgG3

**Clonality:** Monoclonal

**Buffer:** 1.0 mL of PBS containing 0.09% Sodium Azide as preservative.

**Conjugation:** FITC

**Storage:** Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.  
This product is photosensitive and should be protected from light.  
Avoid repeated freezing and thawing.  
Shelf life: one year from despatch.



[View online »](#)