

## Product datasheet for **AM08210RC7-N**

### Mouse IgG3 Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG3
Clonality:	Monoclonal
Buffer:	1.0 mL of PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.
Conjugation:	PE-Cy7
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)