

## Product datasheet for AM08211AC7-N

## 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436

OriGene Technologies, Inc.

Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Mouse IgA Isotype Control**

**Product data:** 

**Product Type:** Isotype Controls

**Applications:** FC

Host: Mouse Isotype: IgA

Clonality: Monoclonal

**Buffer:** 1.0 mL of PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.

**Conjugation:** APC-Cy7

**Storage:** Store the antibody undiluted at 2-8°C.

**DO NOT FREEZE!** 

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing. Shelf life: one year from despatch.