

Product datasheet for **AM08211FC-N**

Mouse IgA Isotype Control

Product data:

Product Type: Isotype Controls

Applications: FC

Host: Mouse

Isotype: IgA

Clonality: Monoclonal

Buffer: 1.0 mL of PBS containing 0.09% Sodium Azide as preservative.

Conjugation: FITC

Storage: Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.
This product is photosensitive and should be protected from light.
Avoid repeated freezing and thawing.
Shelf life: one year from despatch.



[View online »](#)