

Product datasheet for **AM08211RP-N**

Mouse IgA Isotype Control

Product data:

| | |
|---------------|--|
| Product Type: | Isotype Controls |
| Applications: | FC |
| Host: | Mouse |
| Isotype: | IgA |
| Clonality: | Monoclonal |
| Buffer: | 1.0 mL of PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. |
| Conjugation: | PE |
| Storage: | Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch. |



[View online »](#)