

## Product datasheet for **AM08212BT-N**

### Rat IgG1 Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	E, FC
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody is specific to Keyhole limpet hemocyanin (KLH).
Buffer:	PBS containing 0.09% Sodium Azide as preservative.
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)