

Product datasheet for **AM08212PU-N**

Rat IgG1 Isotype Control

Product data:

Product Type: Isotype Controls

Applications: E, FC

Host: Rat

Isotype: IgG1

Clonality: Monoclonal

Specificity: This antibody is specific to Keyhole limpet hemocyanin (KLH).

Buffer: 100 mM Borate buffered saline, pH 8.2.
No preservatives or amine-containing buffer salts added.

Concentration: 0.5 mg/ml

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.
Avoid repeated freezing and thawing.
Shelf life: one year from despatch.



[View online »](#)