

Product datasheet for **AM08213BT-N**

Rat IgG2a Isotype Control

Product data:

Product Type: Isotype Controls

Applications: E, FC

Host: Rat

Isotype: IgG2a

Clonality: Monoclonal

Specificity: This antibody is specific to Keyhole limpet hemocyanin (KLH).

Buffer: PBS containing 0.09% Sodium Azide as preservative.

Concentration: 0.5 mg/ml

Conjugation: Biotin

Storage: Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.
This product is photosensitive and should be protected from light.
Avoid repeated freezing and thawing.
Shelf life: one year from despatch.



[View online »](#)