

Product datasheet for **AM08213FC-N**

Rat IgG2a Isotype Control

Product data:

Product Type:	Isotype Controls
Applications:	E, FC
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody is specific to Keyhole limpet hemocyanin (KLH).
Buffer:	PBS containing 0.09% Sodium Azide as preservative.
Concentration:	0.5 mg/ml
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C. Do NOT freeze. This product is photosensitive and should be protected from light. Shelf life: one year from despatch.



[View online »](#)