

## Product datasheet for **AM08213PU-N**

### Rat IgG2a Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	E, FC
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody is specific to Keyhole limpet hemocyanin (KLH).
Buffer:	100 mM Borate buffered saline, pH 8.2. No preservatives or amine-containing buffer salts added.
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)