

## Product datasheet for **AM08214FC-N**

### Rat IgG2b Isotype Control

#### Product data:

**Product Type:** Isotype Controls

**Applications:** FC

**Host:** Rat

**Isotype:** IgG2b

**Clonality:** Monoclonal

**Specificity:** This antibody is specific to Keyhole limpet hemocyanin (KLH).

**Buffer:** PBS containing 0.09% Sodium Azide as preservative.

**Concentration:** 0.5 mg/ml

**Conjugation:** FITC

**Storage:** Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)