

## Product datasheet for **AM08214PU-N**

### Rat IgG2b Isotype Control

#### Product data:

**Product Type:** Isotype Controls

**Applications:** E, FC

**Host:** Rat

**Isotype:** IgG2b

**Clonality:** Monoclonal

**Specificity:** This antibody is specific to Keyhole limpet hemocyanin (KLH).

**Buffer:** 100 mM Borate buffered saline, pH 8.2.  
No preservatives or amine-containing buffer salts added.

**Concentration:** 0.5 mg/ml

**Conjugation:** Unconjugated

**Storage:** Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.  
Avoid repeated freezing and thawing.  
Shelf life: one year from despatch.



[View online »](#)