

## Product datasheet for **AM08214RC7-N**

### Rat IgG2b Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Rat
Isotype:	IgG2b
Clonality:	Monoclonal
Specificity:	This antibody is specific to Keyhole limpet hemocyanin (KLH).
Buffer:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.
Concentration:	0.1 mg/ml
Conjugation:	PE-Cy7
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)