

## Product datasheet for **AM08215LE-N**

### Rat IgM Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	E, FC, P
Host:	Rat
Isotype:	IgM
Clonality:	Monoclonal
Specificity:	This antibody is specific to NIP.
Buffer:	PBS without preservatives.
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing. Shelf life: one year from despatch.

**Background:** Isotype Control monoclonal antibodies are used to estimate the non-specific binding of target primary antibodies to cell surface antigens. Non-specific binding is due to Fc receptor binding or other protein-protein interactions. Isotype controls show negligible cross-reactivity with cell surface antigens on tissue sections or in cellular preparations. Isotype Controls should be used at identical concentrations and staining conditions as the target primary antibodies.



[View online »](#)