

## Product datasheet for **AM08215RC7-N**

### Rat IgM Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Rat
Isotype:	IgM
Clonality:	Monoclonal
Specificity:	This antibody is specific to NIP.
Buffer:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.
Concentration:	0.1 mg/ml
Conjugation:	PE-Cy7
Storage:	Store the antibody undiluted at 2-8°C.

**DO NOT FREEZE!**

This product is photosensitive and should be protected from light.  
Avoid repeated freezing and thawing.  
Shelf life: one year from despatch.

Background:	Isotype Control monoclonal antibodies are used to estimate the non-specific binding of target primary antibodies to cell surface antigens. Non-specific binding is due to Fc receptor binding or other protein-protein interactions. Isotype controls show negligible cross-reactivity with cell surface antigens on tissue sections or in cellular preparations. Isotype Controls should be used at identical concentrations and staining conditions as the target primary antibodies.
-------------	---



[View online »](#)