

## Product datasheet for **AM09208PU-N**

### HCV NS4 Mouse Monoclonal Antibody [Clone ID: 7E11]

#### Product data:

|                        |                                                                                                                                                                                                                                                                                                                                                                                                       |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:          | Primary Antibodies                                                                                                                                                                                                                                                                                                                                                                                    |
| Clone Name:            | 7E11                                                                                                                                                                                                                                                                                                                                                                                                  |
| Applications:          | ELISA, IHC, WB                                                                                                                                                                                                                                                                                                                                                                                        |
| Recommended Dilution:  | <b>ELISA:</b> React with Human Hepatitis C Virus.<br><b>Western Blot:</b> Use of HCV NS4 antibody (clone 7E11) at 0.5 µg/ml will allow visualization of 0.5 µg/lane of synthetic NS-4 peptide and 0.1 µg/lane recombinant Chimeric HCV polyprotein.<br><b>Immunohistochemistry:</b> Immunoperoxidase-Avidin-Biotin (ABC) assay for Formalin-fixed paraffin embedded tissue section using diluted MAb. |
| Reactivity:            | Human                                                                                                                                                                                                                                                                                                                                                                                                 |
| Host:                  | Mouse                                                                                                                                                                                                                                                                                                                                                                                                 |
| Isotype:               | IgG1                                                                                                                                                                                                                                                                                                                                                                                                  |
| Clonality:             | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                            |
| Immunogen:             | Synthetic polypeptide in which overlapping sequences of at least 90 residues in length were selected from the immunodominant NS-4 regions of Chinese HCV strains.                                                                                                                                                                                                                                     |
| Specificity:           | The HCV NS4 antibody (clone 7E11) reacts with synthetic NS-4a Protein and Recombinant Chimeric HCV Polyprotein.<br>No cross-reaction with HCV Capsid region and other non-structure region.                                                                                                                                                                                                           |
| Formulation:           | 0.01M PBS, pH 7.0 without preservatives.<br>State: Aff - Purified<br>State: Lyophilized purified Ig fraction.                                                                                                                                                                                                                                                                                         |
| Reconstitution Method: | Restore with Double distilled water to adjust the final concentration to 1.0 mg/ml.                                                                                                                                                                                                                                                                                                                   |
| Purification:          | Affinity Chromatography on Protein G.                                                                                                                                                                                                                                                                                                                                                                 |
| Conjugation:           | Unconjugated                                                                                                                                                                                                                                                                                                                                                                                          |
| Storage:               | Store the antibody at -20°C.<br>Avoid repeated freezing and thawing.                                                                                                                                                                                                                                                                                                                                  |
| Stability:             | Shelf life: one year from despatch.                                                                                                                                                                                                                                                                                                                                                                   |



[View online »](#)

### Background:

NS4 A/B are two of the seven nonstructural (NS) proteins making up the HCV polyprotein. The NS proteins recruit the viral genome into an RNA replication complex, which is associated with rearranged cytoplasmic membranes. NS4A acts as a cofactor with the NS3 serine protease and stabilizes its folding. The NS3-NS4A complex is essential for the activation of the latter and allows membrane anchorage of NS3.

HCV is a positive, single-stranded RNA virus in the Flaviviridae family. The genome is approximately 10,000 nucleotides and encodes a single polyprotein of about 3,000 amino acids.

HCV is responsible for a large proportion of worldwide chronic viral hepatitis. Most of these infections develop into chronic hepatitis, which often progresses to liver cirrhosis and hepatocellular carcinoma. At present, (unlike hepatitis A and B), there is no vaccine to prevent hepatitis C infection.

The hepatitis C virus (HCV) nonstructural protein 4B (NS4B) is a relatively hydrophobic 27-kDa protein. The 4A protein has a molecular weight of 6 kDa.

### Synonyms:

Hepatitis C Virus, NS4B, Non-structural protein NS4A