

Product datasheet for **AM09212PU-N**

HCV NS4 Mouse Monoclonal Antibody [Clone ID: 497]

Product data:

Product Type:	Primary Antibodies
Clone Name:	497
Applications:	ELISA, WB
Recommended Dilution:	ELISA: React with Human Hepatitis C Virus. Western Blot: Use of HCV NS4 antibody (clone 497) at 0.5 µg/ml will allow visualization of 0.5 µg/lane of synthetic NS-4 peptide and 0.1 µg/lane recombinant Chimeric HCV polyprotein.
Reactivity:	Human
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Immunogen:	Purified recombinant Chimeric HCV Polyprotein (555 a.a.)
Specificity:	The HCV NS4 antibody (clone 497) reacts with synthetic NS-4a Protein (aa 1689-1735) and Recombinant Chimeric HCV Polyprotein (60 kDa). No cross reaction with Synthetic recombinant Capsid Protein C + Envelope Protein M (core) (aa 1-142), synthetic Capsid Protein C (aa 1-61), recombinant NS-3 Protein (aa 1252-1477) and synthetic NS-3 Protein (aa 1378-1458).
Formulation:	0.01M PBS, pH 7.0 without preservatives. State: Aff - Purified State: Lyophilized purified Ig fraction.
Reconstitution Method:	Restore with Double distilled water to adjust the final concentration to 1.0 mg/ml.
Purification:	Affinity Chromatography on Protein G.
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)

Background:

NS4 A/B are two of the seven nonstructural (NS) proteins making up the HCV polyprotein. The NS proteins recruit the viral genome into an RNA replication complex, which is associated with rearranged cytoplasmic membranes. NS4A acts as a cofactor with the NS3 serine protease and stabilizes its folding. The NS3-NS4A complex is essential for the activation of the latter and allows membrane anchorage of NS3.

HCV is a positive, single-stranded RNA virus in the Flaviviridae family. The genome is approximately 10,000 nucleotides and encodes a single polyprotein of about 3,000 amino acids.

HCV is responsible for a large proportion of worldwide chronic viral hepatitis. Most of these infections develop into chronic hepatitis, which often progresses to liver cirrhosis and hepatocellular carcinoma. At present, (unlike hepatitis A and B), there is no vaccine to prevent hepatitis C infection.

The hepatitis C virus (HCV) nonstructural protein 4B (NS4B) is a relatively hydrophobic 27-kDa protein. The 4A protein has a molecular weight of 6 kDa.

Synonyms:

Hepatitis C Virus, NS4B, Non-structural protein NS4A