

Product datasheet for **AM09388PU-S**

CBR1 Mouse Monoclonal Antibody [Clone ID: AT4E12]

Product data:

Product Type:	Primary Antibodies
Clone Name:	AT4E12
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western blot (1:500 - 1:1000).
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Recombinant human CBR1 (1-277 aa) purified from E. coli
Specificity:	The antibody recognizes human CBR1. Other species not tested.
Formulation:	PBS, pH 7.4 containing 0.02% Sodium Azide and 10% Glycerol State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein-G affinity chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	carbonyl reductase 1
Database Link:	Entrez Gene 873 Human P16152



[View online »](#)

Background:

Carbonyl reductase 1 (CBR1) is a NADPH-dependent, monomeric, and cytosolic enzyme belonging to a family of short-chain dehydrogenases/reductases. This protein consists of 277 amino acid residues and is widely distributed in human tissues such as liver, epidermis, stomach, small intestine, kidney, neuronal cells, and smooth muscle fiber. CBR1 metabolizes many toxic environmental quinones and pharmacological relevant substrates such as the anticancer drug, doxorubicin. The best substrates of CBR1 are quinones, including ubiquinone-1 and tocophrolquinone (vitamin E).

Synonyms:


CBR1, CBR, CRN

Protein Families:


Druggable Genome

Protein Pathways:

Arachidonic acid metabolism, Metabolic pathways

Product images:

The cell lysates (40ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human CBR1 antibody (1:1,000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.



Western blot: The cell lysates (40ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human CBR1 antibody (1:1,000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.