

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for AM09392PU-S

HPRT (HPRT1) Mouse Monoclonal Antibody [Clone ID: AT2G8]

Product data:

Product Type:	Primary Antibodies
Clone Name:	AT2G8
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western blot (1/500-1/5000). Recommended starting dilution is 1/1000.
Reactivity:	Human
Host:	Mouse
lsotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human HPRT (1-218aa) purified from E. coli
Specificity:	The antibody recognizes Human HPRT. Other species not tested.
Formulation:	PBS, pH 7.4 containing 0.02% Sodium Azide and 10% Glycerol State: Purified State: Liquid purified lg fraction
Concentration:	lot specific
Purification:	Protein-G affinity chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	hypoxanthine phosphoribosyltransferase 1
Database Link:	<u>Entrez Gene 3251 Human</u> <u>P00492</u>
Background:	HPRT is a transferase, which catalyzes conversion of hypoxanthine to inosine monophosphate and guanine to guanosine monophosphate via transfer of the 5- phosphoribosyl group from 5-phosphoribosyl 1-pyrophosphate. HPRT, which acts as a catalyst in the reaction between guanine and phosphoribosyl pyrophosphate to form GMP, functions primarily to salvage purines from degraded DNA to renewed purine synthesis.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


	HPRT (HPRT1) Mouse Monoclonal Antibody [Clone ID: AT2G8] – AM09392PU-S	
Synonyms:	HGPRTase, HGPRT	

Druggable Genome, Stem cell - Pluripotency

Protein Families:

Protein Pathways: Drug metabolism - other enzymes, Metabolic pathways, Purine metabolism

Product images:



Western blot analysis: The extracts of 293T (35ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human HPRT (1:500-1:5000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US