

Product datasheet for **AM10167BT-N**

Fluorescein (FITC) Mouse Monoclonal Antibody [Clone ID: 6A4]

Product data:

| | |
|-----------------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | 6A4 |
| Applications: | ELISA, IHC, WB |
| Recommended Dilution: | ELISA: 1/1000-1/10000. Western Blot: 1/500-1/5000 Southern Blot: 1/500-1/5000. <i>In Situ</i> hbridization. Immunohistochemistry on Frozen Sections: 1/200-1/2000. Immunohistochemistry on Paraffin Sections: 1/200-1/2000. |
| Host: | Mouse |
| Isotype: | IgM |
| Clonality: | Monoclonal |
| Immunogen: | Fluorescein (FITC) |
| Specificity: | Reacts with <i>free</i> and <i>bound</i> Fluorescein |
| Formulation: | PBS, pH 7.4 Label: Biotin State: Liquid purified Ig fraction Stabilizer: 1% BSA, 50% Glycerol Preservative: 0.02% Thiomersal |
| Conjugation: | Biotin |
| Storage: | Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |



[View online »](#)