

Product datasheet for **AM10174SU-N**

Ganglioside OAcGD3 Mouse Monoclonal Antibody [Clone ID: 7H2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	7H2
Applications:	ELISA
Recommended Dilution:	ELISA: 1/5-1/50. Thin-layer chromatography (TLC): 1/5-1/20.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG3
Clonality:	Monoclonal
Immunogen:	Ganglioside-OAcGD3
Specificity:	Reacts with OAcGD3 Ganglioside: OAcNeuAc alpha 2-8NeuAc alpha 2-3Gal beta 1-4Glc beta 1-1Cer.
Formulation:	State: Supernatant State: Lyophilized Cell Culture Supernatant Preservative: None
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)