

Product datasheet for **AM20227PU-N**

Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: HP6061]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: HP6061]
Clone Name:	HP6061
Applications:	ELISA
Recommended Dilution:	ELISA: 2.5-10 µg/mL This antibody is best used as a Capture antibody.
Reactivity:	Human
Host:	Mouse
Isotype:	IgM
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.



[View online »](#)