

Product datasheet for **AM20228AP-N**

Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: HP6029]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: HP6029]
Clone Name:	HP6029
Applications:	ELISA
Recommended Dilution:	ELISA: 1/500-1/1,000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Formulation:	50mM Tris/1mM MgCl ₂ /50% Glycerol, pH 8.0 containing 0.09% Sodium Azide as preservative. Label: AP State: Liquid purified Ig fraction. Label: Alkaline Phosphatase
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)