

Product datasheet for **AM20228BT-N**

Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: HP6029]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: HP6029]
Clone Name:	HP6029
Applications:	ELISA
Recommended Dilution:	ELISA: 1/5,000-1/10,000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C.



[View online »](#)