

Product datasheet for **AM20228HR-N**

Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: HP6029]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: HP6029]
Clone Name:	HP6029
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1,000-1/2,000
Reactivity:	Human
Host:	Mouse
Immunogen:	Unknown.
Isotype:	IgG2a
Formulation:	50% Glycerol/50% PBS, pH 7.4 without preservatives. Label: HRP State: Liquid purified Ig fraction. Preservative: No preservatives added. Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)