

Product datasheet for **AM20228PU-N**

Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: HP6029]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: HP6029]
Clone Name:	HP6029
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	ELISA: Purified (UNLB) antibody ≤ 1 ug/mL. Western blotting Dot- and slot-immunoblotting. Immunocytochemistry. Immunohistochemistry on Ffrozen Sections.
Reactivity:	Human
Host:	Mouse
Immunogen:	Unknown.
Isotype:	IgG2a
Formulation:	100 mM Borate buffered saline, pH 8.2 without preservatives or amine-containing buffer salts. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C



[View online »](#)