

## Product datasheet for **AM20229AP-N**

### Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: B3102E8]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: B3102E8]
Clone Name:	B3102E8
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/500-1/1,000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	50mM Tris/1mM MgCl <sub>2</sub> /50% Glycerol, pH 8.0 containing 0.09% Sodium Azide as preservative. Label: AP State: Liquid purified Ig fraction. Label: Alkaline Phosphatase
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)