

Product datasheet for AM20229BT-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: B3102E8]

Product data:

Product Type: Secondary Antibodies

Product Name: Human IgE (Epsilon chain specific) Mouse Monoclonal Antibody [Clone ID: B3102E8]

Clone Name: B3102E8
Applications: ELISA

Recommended Dilution: ELISA: 1/1,000-1/2,000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Formulation: PBS containing 0.09% Sodium Azide as preservative

Label: Biotin

State: Liquid purified Ig fraction.

Concentration: lot specific
Conjugation: Biotin

Storage: Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing.