

Product datasheet for **AM20238PU-S**

Sialylated Mucin-Type-O-Glycans Mouse Monoclonal Antibody [Clone ID: FW17]

Product data:

Product Type:	Primary Antibodies
Clone Name:	FW17
Applications:	ELISA, IF, IHC
Recommended Dilution:	ELISA. Immunocytochemistry. Immunohistochemistry on Paraffin Sections.
Reactivity:	Human
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Immunogen:	Purified mucus glycoproteins (Mucins) from amniotic fluid
Specificity:	Recognizes Sialylated mucin-type O-glycans, oncofetal antigen.
Formulation:	PBS, pH 7.4 State: Purified State: Lyophilized purified Ig fraction
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml
Purification:	Protein G Chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.



[View online »](#)