

Product datasheet for **AM20247AF-N**

Human IgA2 Mouse Monoclonal Antibody [Clone ID: NI 512 (A89-038)]

Product data:

| | |
|------------------------|--|
| Product Type: | Secondary Antibodies |
| Product Name: | Human IgA2 Mouse Monoclonal Antibody [Clone ID: NI 512 (A89-038)] |
| Clone Name: | NI 512 (A89-038) |
| Applications: | ELISA, IF, IP, WB |
| Recommended Dilution: | To identify the presence of IgA2 in Human serum, other body fluids, cell and tissue substrates and to determine its concentration in techniques as Radioimmunoassay, ELISA, indirect Immunoperoxidase and Immunofluorescence staining of cytoplasmic IgA2, Haemagglutination and Immunoblotting. <u>General Recommended Dilutions:</u> Histochemical Use: 1/50-1/200 ELISA: from 1/500 upwards. Western blot: from 1/1000 upwards. |
| Reactivity: | Human |
| Host: | Mouse |
| Immunogen: | Highly purified monoclonal IgA2 isolated from Human serum. |
| Isotype: | IgG1 |
| Formulation: | State: Ascites State: Lyophilized (Delipidated, heat inactivated) stable whole ascites without foreign proteins. |
| Reconstitution Method: | Restore by adding 0.5 ml sterile distilled water. |
| Concentration: | 1.0 mg/ml |
| Conjugation: | Unconjugated |
| Storage: | Store lyophilized at 2-8°C and reconstituted at 2-8°C for one week or (in aliquots) at -20°C for longer. Avoid Repeated thawing and freezing. |



[View online »](#)