

## Product datasheet for **AM20255AF-N**

### Human IgG3 (F(ab)2 specific) Mouse Monoclonal Antibody [Clone ID: NI 330 (HP 6050)]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG3 (F(ab)2 specific) Mouse Monoclonal Antibody [Clone ID: NI 330 (HP 6050)]
Clone Name:	NI 330 (HP 6050)
Applications:	ELISA, IF, IP, WB
Recommended Dilution:	To identify the presence of IgG3 in human serum, other body fluids, cell and tissue substrates and to determine its concentration in techniques as ELISA, indirect immunoperoxidase and indirect immunofluorescence staining of cytoplasmic IgG3, and immunoblotting. <u>General Recommended Dilutions:</u> Histochemical Use: 1/50-1/200. ELISA: from 1/500 upwards. Western blot: from 1/1000 upwards.
Reactivity:	Human
Host:	Mouse
Immunogen:	Highly purified monoclonal IgG1k isolated from pooled human serum
Isotype:	IgG1
Formulation:	State: Ascites State: Lyophilized (Delipidated, heat inactivated) stable ascites without preservatives or foreign proteins
Reconstitution Method:	Restore by adding 0.5 ml sterile distilled water
Concentration:	1.0 mg/ml
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at 2-8°C for one week or (in aliquots) at -20°C for longer. Avoid Repeated thawing and freezing.



[View online »](#)