

Product datasheet for AM20258FC-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Human IgG (Isotype Specific) Mouse Monoclonal Antibody [Clone ID: NI 335 and NI 343]

Product data:

Product Type: Secondary Antibodies

Product Name: Human IgG (Isotype Specific) Mouse Monoclonal Antibody [Clone ID: NI 335 and NI 343]

Clone Name: NI 335 and NI 343 Applications: ELISA, IF, IP, WB

Recommended Dilution: To identify the presence of IgG in human serum, other body fluids, cell and tissue substrates

and to determine its concentration in techniques as ELISA, indirect immunoperoxidase staining of cytoplasmic IgG, and immunoblotting using a peroxidase labelled monoclonal

antibody against FITC.

<u>General Recommended Dilutions:</u> Histochemical Use: 1/10-1/50 ELISA: from 1/200 upwards

Western blot: from 1/400 upwards

Reactivity: Human Host: Mouse

Immunogen: Highly purified monoclonal IgG isolated from human serum.

Isotype: IgG1

Formulation: PBS, pH 7.2 without preservatives and foreign proteins

Label: FITC

State: Lyophilized purified IgG fraction Label: Fluorescein Isothiocyanate Absorption emission: 492 nm / 515 nm Molar radio: Fluorochrome/IgG ~1.1

Reconstitution Method: Restore by adding 0.5 ml sterile distilled water

Concentration: 0.4 mg/ml

Conjugation: FITC

Storage: Store the antibody lyophilized at 2-8°C for one month or (in aliquots) at -20°C for longer.

Diluted ascites should be stored at 2-8°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation.

Avoid repeated thawing and freezing.

