

Product datasheet for **AM20261SU-N**

Human IgA (Secretory component) Mouse Monoclonal Antibody [Clone ID: NI 194-4 (A89-039)]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA (Secretory component) Mouse Monoclonal Antibody [Clone ID: NI 194-4 (A89-039)]
Clone Name:	NI 194-4 (A89-039)
Applications:	ELISA, IF, IHC, IP, WB
Recommended Dilution:	To identify the presence of free or bound Secretory component in Human serum, other body fluids, cell and tissue substrates and to determine its concentration in techniques as Radioimmunoassay, ELISA, Indirect Immunoperoxidase and Immunofluorescence staining, Haemagglutination and Immunoblotting. <i>General Recommended Dilutions:</i> Histochemical Use: from 1/50-1/200. ELISA: form 1/500 upwards. Western blot: from 1/1000 upwards.
Reactivity:	Human
Host:	Mouse
Immunogen:	Highly purified Secretory component isolated from Human milk.
Isotype:	IgG1
Formulation:	State: Ascites State: Lyophilized (Delipidated and heat inactivated) stable whole ascites
Reconstitution Method:	Restore by adding 0.5 ml sterile distilled water. Dilutions may be prepared by adding PBS, pH 7.2
Concentration:	1.0 mg/ml
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)