

## Product datasheet for **AM20437PU-N**


### ART3 (29-139) Mouse Monoclonal Antibody [Clone ID: 3A2]

#### Product data:


Product Type:	Primary Antibodies
Clone Name:	3A2
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry on Paraffin Sections:</b> 5 µg/ml. <b>Western Blot.</b>
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	ART3 (NP_001170, 29 a.a. ~ 139 a.a) partial recombinant protein with GST tag
Specificity:	Recognizes Adp-ribosyltransferase 3 (ART3)
Formulation:	PBS, pH 7.4 State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	ADP-ribosyltransferase 3
Database Link:	<a href="#">Entrez Gene 419 Human Q13508</a>
Synonyms:	FLJ26404; mono-ADP-ribosyltransferase; TMART
Protein Families:	Transmembrane



[View online »](#)

**Product images:**

ELISA: Detection limit for recombinant GST tagged ART3 is approximately 1 ng/ml as a capture antibody.



Human Testis: Formalin-Fixed, Paraffin-Embedded (FFPE)