

Product datasheet for AM20625PU-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Growth Hormone (GH1) Mouse Monoclonal Antibody [Clone ID: GA-8]

Product data:

Product Type: Primary Antibodies

Clone Name: GA-8
Applications: IHC

Recommended Dilution: Immunohistochemistry on paraffin sections: 2 - 4 μg/ml.

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Recombinant hGH.

Specificity: This antibody reacts to Somatotropin.

Formulation: 1.2 % sodium acetate, with 2 mg BSA and 0.01 mg sodium azide as preservative.

State: Purified

State: Lyphilized purified Ig fraction

Reconstitution Method: Restore with 1.2% sodium acetate or neutral PBS

Concentration: 0,1 mg/ml (after reconstitution with PBS)

Purification: Affinity chromatography

Conjugation: Unconjugated

Storage: Prior to reconstitution store at -20°C.

Following reconstitution store undiluted at 2-8°C for one month

or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: growth hormone 1

Database Link: Entrez Gene 2688 Human

P01241





Growth Hormone (GH1) Mouse Monoclonal Antibody [Clone ID: GA-8] - AM20625PU-N

Background: Growth Hormone (GH) is mapped to 17q22-q24. Human growth hormone has a molecular

mass of 22,005 and contains 191 amino acid residues with 2 disulfide bridges. It binds two receptor molecules and thereby induces signal transduction through receptor dimerization. At high concentrations, GH acts as an antagonist because of a large difference in affinities at

the respective binding sites.

Synonyms: Pituitary growth hormone, Growth hormone 1, HGH

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway, Neuroactive ligand-

receptor interaction