

## Product datasheet for **AM20654PU-N**

### CDK1 Mouse Monoclonal Antibody [Clone ID: IMD-34]

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | Primary Antibodies   |
| Clone Name:       | IMD-34   |
| Applications:     | IHC, WB  |
| Reactivity:       | Chicken, Human, Mouse, Rabbit, Rat   |
| Host:             | Mouse  |
| Isotype:          | IgG2a  |
| Clonality:        | Monoclonal   |
| Formulation:      | State: Purified  |
| Conjugation:      | Unconjugated   |
| Gene Name:        | cyclin-dependent kinase 1  |
| Database Link:    | <a href="#">Entrez Gene 12534 Mouse</a> <a href="#">Entrez Gene 54237 Rat</a> <a href="#">Entrez Gene 983 Human P06493</a> |
| Synonyms:         | CDK1, CDC2, CDC28A, CDKN1, P34CDC2, p34 protein kinase   |
| Protein Families: | Druggable Genome, Protein Kinase, Stem cell - Pluripotency   |
| Protein Pathways: | Cell cycle, Gap junction, Oocyte meiosis, p53 signaling pathway, Progesterone-mediated oocyte maturation                   |



[View online »](#)