

Product datasheet for **AM20780AG-L**

DYKDDDDK Epitope Tag Mouse Monoclonal Antibody [Clone ID: 3B9]

Product data:

Product Type:	Primary Antibodies
Clone Name:	3B9
Applications:	IP
Recommended Dilution:	Immunoprecipitation.
Reactivity:	All Species
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Synthetic peptide containing epitope DYKDDDDK (KLH-coupled).
Specificity:	This antibody reacts to DYKDDDDK Epitope Tag.
Formulation:	50% slurry in storage buffer (1× PBS, pH 7.4, containing 0.09% sodium azide). <u>Recommended elution buffer:</u> 0.2 M Glycine, pH 2.5 Label: Agarose State: Liquid purified Ig fraction
Conjugation:	Agarose
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.
Background:	Anti-DYKDDDDK-Tag Mouse mAb (Agarose Conjugated) is a monoclonal anti-DYKDDDDK antibody covalently linked to agarose; the agarose enables immunoprecipitation (IP) of DYKDDDDK tagged proteins or co-immunoprecipitation (Co-IP) of their interacting partners.
Synonyms:	FLAG Epitope Tag, ECS Epitope Tag, FLAG-tag, ECS-tag, D-tag



[View online »](#)

- Note:** Protocol: Immunoprecipitation procedure
The work can be performed in 1.5 ml micro-centrifuge tubes or in spin columns.
1. Thoroughly resuspend the Anti-DYKDDDDK Agarose by inverting the tube or vial several times.
 2. Add 20-50 μ l 50% slurry of Anti-DYKDDDDK Agarose into cell lysate using a wide-bore pipette tip. Note: the lysate should be fresh, and for a well expressed tagged protein, 200 μ l lysate (200-500 g total protein) usually yields a good IP result.
 3. Incubate with gentle mixing for 2 h to overnight at 4°C.
 4. Wash the beads with 1 ml TBS buffer or lysis buffer, such as RIPA (50 mM Tris HCl, pH 7.4, 150 mM NaCl, 1 mM EDTA, 1% NP-40, 0.5% sodium deoxycholate), centrifuge for 3 min at 2,000 \times g, and discard the supernatant. Wash 3 times, avoid losing beads during washes.
 5. Elution of the DYKDDDDK tagged protein.
Option 1. Elution with elution buffer.
Add 30-50 μ l elution buffer to the beads, gently tap the tube to mix well, immediately centrifuge for 3 min, transfer the supernatant very carefully to a fresh tube (Avoid transferring any beads).
Note: Neutralize the eluant immediately by add 1 μ l of 1.5 M Tris, pH 9.0 per 20 μ l Elution buffer.

Option 2. Elution with DYKDDDDK peptide
Add 30-50 μ l DYKDDDDK peptide solution (100 g/ml DYKDDDDK peptide in TBS buffer), gently tap the tube to mix well, incubate for 10 min, centrifuge for 3 min, and transfer the supernatant to a fresh tube. TBS buffer: 50 mM Tris HCl, 150 mM NaCl, pH 7.4.

Option 3. Elution with SDS loading buffer
Add 30 μ l 2 SDS loading buffer, gently tap the tube to mix well, boil at 100°C for 5 min, centrifuge for 3 min, transfer the supernatant to a fresh tube.
Note: in this case, the supernatant contains not only the binding proteins, but also IgG (heavy and light chains).
 6. Prepare SDS-PAGE gel for western blotting or proceed to other assays