

## Product datasheet for **AM20787PU-L**

### **GST-Tag Rat Monoclonal Antibody [Clone ID: KT4]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	KT4
<b>Applications:</b>	ELISA, WB
<b>Recommended Dilution:</b>	Suitable for ELISA (1 µg/ml) and Western Blotting (1 µg/ml).
<b>Host:</b>	Rat
<b>Isotype:</b>	IgG2a
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	A fusion protein containing GST
<b>Specificity:</b>	Recognizes Glutathione-S-transferase (GST).
<b>Formulation:</b>	PBS containing 0.09% Sodium Azide as a preservative. State: Aff - Purified State: Liquid purified Ig fraction.
<b>Concentration:</b>	lot specific
<b>Purification:</b>	Protein G Chromatography.
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.
<b>Database Link:</b>	<a href="#">P08515</a>
<b>Synonyms:</b>	Glutathione S-transferase Tag, GST26-Tag



[View online »](#)

**Product images:**

