

Product datasheet for **AM20787PU-S**

GST-Tag Rat Monoclonal Antibody [Clone ID: KT4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT4
Applications:	ELISA, WB
Recommended Dilution:	Suitable for ELISA (1 µg/ml) and Western Blotting (1 µg/ml).
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	A fusion protein containing GST
Specificity:	Recognizes Glutathione-S-transferase (GST).
Formulation:	PBS containing 0.09% Sodium Azide as a preservative. State: Aff - Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Purification:	Protein G Chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	P08515
Synonyms:	Glutathione S-transferase Tag, GST26-Tag



[View online »](#)

Product images:

