

Product datasheet for **AM20790PU-L**

CD89 (FCAR) Mouse Monoclonal Antibody [Clone ID: MIP68b]

Product data:

Product Type:	Primary Antibodies
Clone Name:	MIP68b
Applications:	WB
Recommended Dilution:	Western Blot: 1 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Soluble CD89
Specificity:	This antibody reacts to CD89.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	Fc fragment of IgA receptor
Database Link:	Entrez Gene 2204 Human P24071
Synonyms:	IgA Fc receptor
Protein Families:	Transmembrane



[View online »](#)