

Product datasheet for **AM20791PU-N**

Lyl1 Rat Monoclonal Antibody [Clone ID: KT43]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT43
Applications:	IP, WB
Recommended Dilution:	Western blotting: 1 µg/ml. Immunoprecipitation: 2 µg/ml.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	LYL1 peptide coupled to carrier protein
Specificity:	This antibody reacts to LYL1.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	lymphoblastomic leukemia 1
Database Link:	Entrez Gene 17095 Mouse P27792
Synonyms:	Protein lyl-1, BHLHA18



[View online »](#)