

Product datasheet for **AM20792PU-N**

Cbfa2t2 Rat Monoclonal Antibody [Clone ID: KT42]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT42
Applications:	IP, WB
Recommended Dilution:	Western blotting: 1 µg/ml. Immunoprecipitation: 2 mg/ml.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2c
Clonality:	Monoclonal
Immunogen:	MTGR1 peptide couple to carrier protein
Specificity:	This antibody reacts to CBFA2T2.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	core-binding factor, runt domain, alpha subunit 2, translocated to, 2 (human)
Database Link:	Entrez Gene 12396 Mouse O70374
Synonyms:	EHT, MTGR1, MTG8-like protein, MTG8-related protein 1



[View online »](#)