

## Product datasheet for **AM20792PU-S**

### **Cbfa2t2 Rat Monoclonal Antibody [Clone ID: KT42]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	KT42
<b>Applications:</b>	IP, WB
<b>Recommended Dilution:</b>	Western blotting: 1 µg/ml. Immunoprecipitation: 2 mg/ml.
<b>Reactivity:</b>	Mouse
<b>Host:</b>	Rat
<b>Isotype:</b>	IgG2c
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	MTGR1 peptide couple to carrier protein
<b>Specificity:</b>	This antibody reacts to CBFA2T2.
<b>Formulation:</b>	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
<b>Concentration:</b>	lot specific
<b>Purification:</b>	Protein A chromatography
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.
<b>Gene Name:</b>	core-binding factor, runt domain, alpha subunit 2, translocated to, 2 (human)
<b>Database Link:</b>	<a href="#">Entrez Gene 12396 Mouse</a> <a href="#">O70374</a>
<b>Synonyms:</b>	EHT, MTGR1, MTG8-like protein, MTG8-related protein 1



[View online »](#)