

Product datasheet for **AM20797PU-L**

Nipbl (Short form) Rat Monoclonal Antibody [Clone ID: KT55]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT55
Applications:	IF, IHC, IP, WB
Recommended Dilution:	Western blot (1 µg/ml). Immunofluorescent staining (10 µg/ml). Immunoprecipitation. Immunohistochemistry (10 µg/ml).
Reactivity:	Human, Mouse
Host:	Rat
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Delangin peptide coupled to carrier protein (BSA)
Specificity:	This antibody reacts to NIPBL (Short form).
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	Nipped-B homolog (Drosophila)
Database Link:	Entrez Gene 71175 Mouse Q6KCD5
Synonyms:	IDN3, Delangin, SCC2 homolo



[View online »](#)